


inspired™ O2FLO

breathing life



High Flow Heated Respiratory Humidifier

User Manual



Vincent Medical Manufacturing Co., Ltd.
Flat/RM B2, 7/F., Hang Fung Industrial
Building, Phase 2, 2G Hok Yuen Street,
Hung Hom, Kowloon, Hong Kong
Tel: +852 2365 5688
Fax: +852 2765 8428



Website: www.inspired-medical.com

Sales Inquiry:

+852 2365 5688

info@inspired-medical.com

Service Support:

+852 2186 1010

service@inspired-medical.com

Vincent Medical Manufacturing Co., Ltd

Thank you for purchasing the inspired™ O2FLO High Flow Heated Respiratory Humidifier (VUN-001).

The inspired™ O2FLO is a medical device. Please follow the instructions for use and intended use in this manual to ensure the device is used properly and in a safe manner.

Please keep this manual for future reference.

This manual and its corresponding intellectual property rights belong to Vincent Medical Manufacturing Co., Ltd. Without written consent, any individual or organization must not copy, modify, or translate any part of this manual.

Contents

1. Warnings and Cautions	1
1.1 Intended Use	7
1.2 Product Features	7
2. inspired™ O2FLO Parts, Features and Specification ...	8
2.1 inspired™ O2FLO Main Unit	8
2.2 Package Content	10
2.3 Optional Medical Accessories (Sold Separately)	10
2.4 Control Panel	13
2.5 Specification	15
3. inspired™ O2FLO Operating Instruction	16
3.1 Preparation before Use	16
3.2 Device Setup	19
3.3 Humidification Treatment	24
3.4 Settings Menu	24
3.5 Alarms	26
3.6 After Use	30
3.7 Data Record Software	31

4. Product Maintenance	35
4.1 Filter Replacement	35
4.2 Main Unit Cleaning	36
4.3 Temperature Probe Cleaning and Disinfection	36
4.4 Troubleshooting	37
4.5 Maintenance Statement	38
4.6 Battery Information	38
4.7 Main Unit Leakage Test	39
5. EMC Information	40
5.1 Declaration - Electromagnetic Emissions	40
5.2 Guidance & Declaration - Electromagnetic Immunity	41
5.3 Guidance & Declaration -- Electromagnetic Immunity Concerning Conducted RF & Radiated RF	42
5.4 Recommended Separation Distances Between Portable and Mobile RF Communications Equipment and the Model VUN-001	43
6. Use and Warranty Period	44



1. Warnings and Cautions



Read this manual before use. Keep this instruction manual available for future reference.

Please read the following carefully before use.

The warning icons detailed below are intended to help the user to operate the device in a safe and correct manner and to prevent injuries to the user and others.

The icons and their meanings are as follows.

Examples of Icons		
	Warning, Electricity	Icon indicates caution (including warning and danger).
	General Prohibition	Icon indicates prohibited actions (what you cannot do).
	General Mandatory	Icon indicates something that is compulsory (must be observed at any time).
	Warning; Flammable Material	Icon indicates caution (including warning and danger).



Warnings



- The device is for use with spontaneously breathing patients and not for life support.
- Continuous monitoring of blood oxygen concentration of the patient is required while using this device.
- The device should only be operated by trained personnel under guidance of a physician.
- The device should be located in a position where ventilation around the unit is not restricted.
- This device can be used in desktop, pole-mounted and wall-mounted configurations. Please ensure that the device is not tilted at an angle greater than 10° from the horizontal plane. Failure to do so may cause the device to overbalance and/or may lead to water spilling out from the chamber.
- This device must be installed with an air inlet filter during use.
- The air inlet filter should be replaced at least once a month in single-patient use, and after every patient. Only approved filters supplied by Vincent Medical should be used.
- When the device is not in use, please place it in its original packaging to keep the air inlet and outlet clean.
- The device should be maintained and serviced only by qualified or trained personnel.
- When the device is used with a supply of oxygen, please pay attention to the following:
 - > Inspect the oxygen tube and ensure it is free of kinks.
 - > Keep device away from sources of ignition and/or open flames. Do not smoke near the device.
 - > Oxygen concentration and flow of delivered gas should only be determined by a physician or trained medical professional.
 - > Oxygen supply shall be standard dry and clean compressed oxygen for medical applications. The minimum oxygen concentration required is 99.5% .
 - > Adjust oxygen flow meter slowly.
 - > Do not turn up the oxygen supply until the device is running.
 - > Shut off the oxygen supply before turning off the device.
- Oxygen leakage may cause serious accidents.
 - > If an oxygen leak is discovered, shut off the supply of oxygen and contact authorized personnel for repair.

- The operator or responsible department should contact Vincent Medical or a Vincent Medical representative:
 - > for assistance, if needed, in setting up, using or maintaining the device, or;
 - > to report unexpected operation or events.
- After using the device, store it away from pets and out of the reach of children. Keep the storage environment clean and free from dust.
- Any liquid spilt or entering into the device may cause damage to the device or injury to patient and/or end user.
- Ensure that the chamber is not tilted in any orientation to avoid water entering the breathing circuit.
- Use sterile water bags for inhalation only.
- The maximum surface temperature of the following parts (at ambient temperature of 28°C) are:
 - Heated Breathing Circuit: 43°C;
 - Applied Part (Nasal Cannula): 42.9°C;
 - Stand-by Button: 59.8°C;
 - Control Knob: 44.4°C;
 - LCD Screen: 42.1°C.



- Do not use the device in the presence of a flammable anaesthetic mixture with air, oxygen or nitrous oxide. Doing so could lead to fire or explosion.
- Do not use the device for any purpose other than what is specified in the manual. Doing so may cause damage to the device or injury to patient and/or end user.
- Do not block or obstruct the air and oxygen inlets or the delivery gas outlet. Blockage or obstruction may lead to device malfunction and/or injury to patient and/or end user.
- Do not use any parts or accessories other than those listed in this manual.
- Do not place liquid containers on or near the top of the device.
- Do not attempt to disassemble, repair or modify the unit. Doing so may cause damage to the device or injury to patient and/or end user.
- To avoid contamination of the device, do not blow into the air inlet or outlet.
- Do not fill water chamber beyond the maximum water level indicated on the chamber. Doing so may cause water to spill into the breathing circuit. If water level exceeds the maximum, please replace the chamber with a new one.
- Do not fill the humidification chamber manually from the circuit port.
- Do not pour water with a temperature of over 37°C or below 10°C into the chamber.
- Do not kink or block the water supply tube.




























- Do not apply excessive force to the device.
- Do not drop or subject the device to strong impacts or shocks.
- Avoid exposing the device to direct sunlight and/or high temperatures.
- Do not allow children or infants to tamper with or operate the device.
- Do not immerse the machine in liquid/water.
- Do not use corrosive solvents to clean the surface of the device.
- Do not use the device if it is not working properly. If an abnormality in device function is detected, stop using the device and have a qualified technician or engineer inspect the device before attempting to use it again.
- Do not attempt to service the device while it is in use. Doing so may cause damage to the device or injury to patient and/or end user.
- To prevent burns, do not touch the heater plate when the device is running and do not touch for 30min after running.
- To prevent burns, do not continuously touch the bottom of the chamber for more than 1 minute when it is hot.
- Do not use the device when room temperature is below 18°C or above 28°C.



- Do not use the power cord or the power plug if it is damaged, and do not plug the power cord into a loose electrical outlet/socket. Doing so may lead to electric shock.
- Keep power cord away from heated or hot surfaces as this may damage the power cord. Use of a damaged power cord may lead to electric shock or serious device malfunction.
- Do not plug the power plug into the electrical outlet with wet hands. Doing so may lead to electric shock.
- Do not wash the machine with water, or splash water to the power source. Doing so may lead to electric shock or serious device malfunction.
- If the machine is not being used for a long period of time, it is recommended that the device is unplugged/ disconnected from the electrical outlet.
- To unplug the power plug from the electrical outlet, do not drag the power cord. Unplug by pulling the power plug. Dragging the power cord may lead to electric shock or serious device malfunction.
- Be sure to turn off the power and unplug the power plug from the electrical outlet before installing, moving, or servicing the device. Failure to do so may lead to electric shock or serious device malfunction.
- Do not excessively stretch, bend or apply excessive force to cables and hoses.

NOTE: The device should be maintained and serviced only by qualified or trained personnel.

Explanation of Symbols

	Refer to instruction manual		Consult instructions for use
	Caution		Non-sterile
	Do not re-use		Do not use if package is damaged
	Date of manufacture		Manufacturer
	Use by date		Serial No.
	Batch code		Catalogue number
	Caution, hot surfaces		EU representative
IP21	Particulate matter and drip proof protection rating IP21		Alternating current
	Humidity limitation		Temperature limit
	Separate collection for EEE (EU only)		CE marking (EU only)
	Type BF Applied Part		Class II equipment
	Atmospheric pressure limitation		Do not roll
	This way up		Fragile, handle with care
	Keep dry		Keep away from sunlight

Safety Classification

- Type of protection against electric shocks: Class II .
- Degree of protection against electric shock: Type BF applied parts (Patients will come in direct contact with nasal cannula during treatment.)
- Degree of protection against harmful ingress of particulate matter and water: IP21 .
- Degree of safety of application in the presence of a flammable anaesthetic mixture with air or oxygen or nitrous oxide: Equipment is not suitable for use in the presence of a flammable anaesthetic mixture with air or oxygen or nitrous oxide.
- Mode of operation: Continuous operation.

Environmental Conditions

This device should not be exposed to excessive vibration, dust, corrosive or explosive gases. During device operation it should be placed in a horizontal position. Suitable environmental conditions for operation are:

- > Ambient temperature: 18 ~ 28°C;
- > Relative humidity: 15% ~ 93% RH, non-condensing;
- > Atmospheric pressure: 86 ~ 106kPa;
- > Altitude: up to 2000m above sea level;
- > Storage and transportation conditions: -20°C ~ 50°C, 15% ~ 93% RH, non-condensing, 86 ~ 106kPa.

Environmental Protection



At the end of service life, the device and its accessories should be disposed of and decommissioned in accordance with local laws and regulations. Contact the appropriate local authorities to determine the proper method of disposal of potentially bio-hazardous parts and accessories.

Overview

The inspired™ O2FLO is a humidifier that delivers warmed and humidified respiratory gases to the patient via a flow generator with a flow rate from 2LPM to 60LPM. The gas delivered is composed of an air-oxygen mixture and is capable of delivering oxygen concentration from 21% to 95%.

1.1 Intended Use

The inspired™ O2FLO is intended to deliver warmed and humidified respiratory gases to spontaneously breathing patients with a flow rate from 2LPM to 60LPM and is capable of providing high flow therapy. The device can be used in the emergency room, outpatient department, inpatient department and other diagnostic and treatment rooms. This device can also be used for homecare.

1.2 Product Features

a. Dew Point Temperature Control

This device allows the operator to set the dew point temperature from 31°C to 37°C and is equipped with highly sensitive temperature probes to monitor patient end temperature. The device also has a number of temperature protection settings to ensure patient safety.

b. Oxygen Concentration Control

This device can be used with an additional oxygen supply. The device accurately displays real-time oxygen concentration of gas mixture delivered. The user can control the oxygen concentration by adjusting the external oxygen valve. The upper and lower limits of oxygen concentration alarm is settable to ensure patient safety.

c. Flow Rate Control

This device allows the operator to set the flow rate from 2LPM to 60LPM. The device monitors the flow rate with a built-in highly sensitive flow rate sensor.

d. Mode Control

The device has two operating modes, namely Adult Mode and Paediatric Mode. The device will identify if the correct breathing circuit is being used (either Adult or Paediatric) to avoid misuse.

e. Display and Control Interface

The device is equipped with a LCD display monitor and control knob to make the device easier to control, and to provide clear display of functions and data.

f. Alarms


The device has an auditory alarm function to alert patient or end user when an alarm condition occurs. The auditory alarm is accompanied by indicator light and information displayed on the monitor.



2. inspired™ O2FLO Parts, Features and Specification



2.1 inspired™ O2FLO Main Unit



















No.	Name	Function
①	Stand-by Button	Press and hold for 2 seconds to turn on or off the inspired™ O2FLO.
②	Control Knob	a. Rotate clockwise or counter-clockwise to highlight settings. b. Press the knob to select the highlighted settings c. When alarm is active, press the knob to mute the alarm audio for 120 seconds. Press again to recover alarm audio.
③	Start/ Pause Button	Press Start/ Pause button to Start/ Pause the running.
④	LCD Display	2.8-inch LCD display. Automatically enters the screen saver mode if there is no operation for 10 min.
⑤	Status LED	a. Green (flashing): Device is connected to mains power, but is not turned on. b. Green (continuous): Normal operation. c. Yellow (flashing): Medium Priority Alarm. d. Red (flashing): High Priority Alarm.
⑥	Power Cord	Connect to mains power supply.
⑦	Bracket	Couple with the Bracket Base, to mount the device onto the pole or wall.
⑧	Heater Plate	Heat the water in the chamber.
⑨	Chamber Lock	Fit the chamber to the right position.
⑩	Annular Tube	Air-oxygen mixture outlet, connect to the chamber inlet.
⑪	Water Level Sensor and Chamber Sensor	Detect water level and presence of chamber, to provide water level alarm and chamber uninstalled alarm.
⑫	Temperature Probe Port (Blue)	Connect to Temperature Probe.
⑬	Heater Connection Cable Port (Red)	Connect to Heater Connection Cable.
⑭	SD Card Slot	Connect to SD card.
⑮	Oxygen Inlet Port	Connect to oxygen supply with oxygen tube.
⑯	Filter Cover	Secure the filter in place.
⑰	Filter	Filter the air taken in.
⑱	Heater Connection Cable	For connection between device and 510-090/510-091 Heated Breathing Circuit.
⑲	Temperature Probe	To detect the chamber outlet temperature, patient end temperature and ambient temperature when used 510-090/510-091 Heated Breathing Circuit.
⑳	Connection Cable	For connection between device and 510-049/510-050 Integrated Heated Breathing Circuit.
㉑	Bracket Base	To be screwed onto pole or wall where the inspired™ O2FLO is to be mounted.





2.2 Package Content

Item	Unit
inspired™ O2FLO	1 PC
Heater Connection Cable	1 PC
Temperature Probe (VUN-200)	1 PC
Bracket Base	1 PC
Filter	5x5 PCS
Oxygen Tube	1 PC
User Manual	1 PC
Quick Guide	1 PC
VUN-300 Connection Cable	1 PC

2.3 Optional Medical Accessories (Sold Separately)

Drawing & Product Name	Product Name	P/N	HC/HU	Maximum Period of Use	Single Use
	VHC60 Autofeed Humidification Chamber	51006488 (JP) 51006136 (US/EU/UK)	HC & HU	14 days	Yes
	Oxygen Tube	51007929 (JP) 51006741 (US/EU/UK)	HC & HU	N/A	No
	Filter	51006764	HC & HU	1 month	No
	Heater Connection Cable	51006736	HC & HU	N/A	No
	VUN-200 Temperature Probe (Used for 510-090 and 510-091)	51006737	HC & HU	500 times	No
	510-090 Adult Heated Breathing Circuit	51006729	HC & HU	14 days	Yes

Drawing & Product Name	P/N	HC/HU	Maximum Period of Use	Single Use
 510-091 Paediatric Heated Breathing Circuit	51006728	HC & HU	14 days	Yes
 510-049 Adult Integrated Heated Breathing Circuit	51006765	HC & HU	14 days	Yes
 510-050 Paediatric Integrated Heated Breathing Circuit	51006766	HC & HU	14 days	Yes
 VUN-300 Connection Cable (Used for 510-049 and 510-050)	51006767	HC & HU	500 times	No
 Nasal Cannula (Adult) VANC-01(S) VANC-02 (M) VANC-03 (L)	51006704 (JP) 51006705 (JP) 51006706 (JP) 51006177 (EU) 51006178 (EU) 51006179 (EU)	HC & HU	14 days for HU 30 days for HC	Yes
 Nasal Cannula (Infant) VINC-01 (S)  VINC-02 (M)  VINC-03 (L)  VINC-04 (XL) 	51005447 (EU) 51005243 (EU) 51005244 (EU) 51005245 (EU)	HU	7 days	Yes

Drawing & Product Name	P/N	HC/HU	Maximum Period of Use	Single Use
 Heated Breathing Circuit Kit (3PIN Adult) 511-090-C60	51007895	HC & HU	14 days	Yes
 Heated Breathing Circuit Kit (3PIN Pediatric) 511-091-C60	51007896	HC & HU	14 days	Yes
 Heated Breathing Circuit Kit (6PIN Adult) 511-049-C60	51007893	HC & HU	14 days	Yes
 Heated Breathing Circuit Kit (6PIN Pediatric) 511-050-C60	51007894	HC & HU	14 days	Yes
Remark: Home Care (HC), Hospital Use (HU).				

2.4 Control Panel



Status LED	<ul style="list-style-type: none"> a. Green (flashing): Device is connected to mains power, but is not turned on. b. Green (continuous): Normal operation. c. Yellow (flashing): Medium Priority alarm. d. Red (flashing): High Priority alarm.
Stand-by Button	Press and hold for 2 seconds to turn the device on or off.
Control Knob	<ul style="list-style-type: none"> a. Rotate clockwise or counter-clockwise to highlight settings. b. Press the control knob to select highlighted setting. c. When alarm is active, press the control knob to mute the alarm audio for 120 seconds. Press again to recover alarm audio.
Start/ Pause Button	Press Start/ Pause button to start the device or to pause it while running.
LCD Display	Screen saver will be displayed after 10 minute if no user operation is detected.

Icons

Icon	Name	Function
	Adult Mode	Select to run in Adult Mode. Adult breathing circuit and adult high flow nasal cannula must be used in Adult Mode.
	Paediatric Mode	Select to run in Paediatric Mode. Paediatric breathing circuit and paediatric high flow nasal cannula must be used in Paediatric Mode.
	Oxygen Concentration	Display oxygen concentration and select to set oxygen concentration alarm limits.
	Dew Point Temperature	Select to set Dew Point Temperature at patient end.
	Flow Rate	Select to set Flow Rate.
	Alarm	When the alarm audio is muted, icon will be displayed as
	SD Card	Icon indicates SD card is inserted. When SD card is not inserted, icon will be displayed as
	Setting	Select to set language, time, alarm volume and to display device information and alarm history.
	Back	Select to return to previous page.
	Lock/Unlock	Select to Lock/Unlock settings change.

2.5 Specification

Power Input	100-240V ~ , 50/60Hz
Power	300VA
Flow Rate Range (Adult Mode)	10~60L/min
Flow Rate Range (Paediatric Mode)	2~25L/min
Oxygen Concentration Accuracy	±5%
Working Pressure	30~80 psi (0.2Mpa~0.55Mpa)
Dew Point Temperature Range (Adult Mode)	31~37°C, ±2°C
Dew Point Temperature Range (Paediatric Mode)	31~34°C, ±2°C
Humidification Output	>33mg/L @ 37°C ; >12mg/L @ 31~36°C
Sound Pressure Level	≤50 dB(A)
Dimensions	243mm x 220mm x 170mm
Net Weight	2.8±0.2kg
Expected Service Life of the Unit	5 years




3. inspired™ O2FLO Operating Instruction



3.1 Preparation before Use

1) Install the autofeed humidification chamber by pressing down the chamber lock and sliding the chamber into position sitting on the heater plate. Once the chamber is in the correct position, the chamber lock will click into place with the chamber inlet connecting directly into the annular tube of the machine.




2) Connect water bag

Push the spike at the end of water supply tube into the port at the bottom of the bag. Open the vent cap on the side of the spike. Hang the water bag above the unit. The chamber will be automatically filled to the required level and maintain at that level.

Note: The unit will beep and the status LED will flash in yellow if there is insufficient water in the chamber. In this situation, check and make sure that the water supply tube is free of kinks and that there is sufficient water in the water bag. Replace as necessary if a kink is found or the water bag is empty.

VHC60: Flow setting vs Typical usage time (2 liter sterile water bag).

L/min	2	5	10	15	20	25	30	35	40	45	50	55	60
hrs	340	136	68	46	34	27	22	19	17	15	14	12	11




- 3) Connect Heated Breathing Circuit to the chamber.
- 4) Connect the heater connection cable and Temperature probe as indicated.
 - Blue Port – Temperature Probe
 - Red Port – Heater connection cable

NOTE:

- Do not allow the heated breathing circuit to come into direct contact with the patient's skin for prolonged periods of time.
- Do not use any accessory which is not recommended by Vincent Medical in this manual.

If an integrated heated circuit with integrated temperature probe is used, VUN-300 Connection Cable shall be used in place of Heater Connection Cable and VUN-200 Temperature Probe. The connection is shown as below:



Warning:

- Replace chamber if water fill exceeds maximum level.
- Disconnect water supply when machine is not running.
- Do not use tap water, saline water or glucose water. Use sterile water for inhalation only.
- Do not fill water through chamber outlet.
- Do not tilt machine.
- Do not use chamber for more than 14 days.
- The chamber are all single use products. Do not reuse or use in multiple patients.


5) Select patient interface

The following table shows the dew point temperature and flow settings associated with the relevant patient interfaces. Select the correct patient interfaces as per the guidance of physician or healthcare professional.


Patient Interface	Model	Dew Point Temperature (°C)	Flow Rate (L/min)
Nasal Cannula (Infant)	VINC-01 (S)	31~34	2~8
	VINC-02 (M)	31~34	2~8
	VINC-03 (L)	31~34	2~10
	VINC-04 (XL)	31~34	2~10
Nasal Cannula (Adult)	VANC-01 (S)	31~37	10~50
	VANC-02 (M)	31~37	10~60
	VANC-03 (L)	31~37	10~60


* In low ambient temperature conditions, the dew point temperature may not reach 37°C when the flow is high. User is advised to reduce the target flow setting to achieve the desired dew point temperature.

Warning: The patient interfaces are all single use products. Do not reuse or use in multiple patients as this may lead to cross-contamination. The patient interfaces must be changed frequently to avoid infection.




3.2 Device Setup




- 1) Plug in the power cord.
- 2) Press and hold the power button  for 2 seconds to start the inspired™ O2FLO. The screen will display the Logo and device information (device name, model, software version, serial number, date and time), and the system will start self-test.

During self-test, check that the following functions are working properly:

- The status LED lights up in red, green and yellow sequentially.
- The screen displays the Logo and device information.
- The device generates a beep.





- 3) The time setting page will appear when it is the first time the inspired™ O2FLO starts or after the replacement of internal battery. After setting the date and time, rotate the control knob to select "Save" and press control knob to continue.





- 4) The screen will show the main control interface. The default operating mode is Adult Mode.



5) Parameter Setting

Step 1: Rotate the control knob to highlight the lock icon .



Press control knob to unlock .



Step 2: Rotate the control knob to highlight the mode icon .


Press the control knob. The icon will appear as yellow .


Rotate the control knob to select Paediatric Mode , then press the control knob to confirm. The icon will then appear as green .


If the control knob is not pressed within 10 seconds to confirm the change, the system will return to its previous settings.

Step 3: Use the control knob in a similar manner to set the Dew Point Temperature  (31- 37°C), the target Flow  (2-60 L/min) and the Oxygen Concentration alarm upper limit O₂ 21⁹⁵/₁₈% (default 95%) and lower limit O₂ 21⁹⁵/₁₈% (default 18%).




- 6) Press the Start/Pause Button  to start running. The real-time dew point temperature, flow rate and oxygen concentration will be shown on the screen. The 4 icons of settings will become blue in running.



- 7) To stop running, press the Start/Pause  Button. If oxygen concentration is lower than 25%, this window will appear, select "Yes" and press the control knob to stop, or rotate the control knob to select "No" and press the control knob to continue.

* When the unit stops running, the run time will stop and reset to zero.




If oxygen supply has not been turned off and oxygen concentration is higher than 25%, a window will appear. User must discontinue oxygen supply before stopping the device from running.



8) Adjusting oxygen concentration levels

Connect the oxygen tube between the flowmeter (connected to oxygen supply) and the oxygen inlet port at the right side of the unit. Push the oxygen tube firmly into this port to ensure a secure connection.

** The flowmeter shall be capable of delivering oxygen flow up to 60LPM, otherwise the oxygen concentration under high flow setting would be limited to low levels.*



Adjust the oxygen flow while monitoring the oxygen concentration O_2 shown on the screen. Refer to the table next page to achieve the target O_2 %.




Flow to Patient (L/min) / Flow-meter (L/min) / O ₂ %	O ₂ %													
	4	6	8	10	15	20	25	30	35	40	45	50	55	60
2	60	47	41	37	32	29	27	26	26	25	25	24	24	24
4		73	60	52	42	37	34	30	30	29	28	27	27	26
6			80	68	52	45	40	37	34	33	31	30	30	29
8				84	63	52	46	42	39	37	35	34	32	31
10					73	61	53	47	44	41	39	37	35	34
15						80	68	61	55	50	47	45	43	41
20							84	74	66	60	56	53	50	47
30									89	80	74	68	64	61
40											91	84	78	74
50													93	87

Attention: When oxygen supply is not connected to the device, please keep the cap of the oxygen inlet port closed.

Warning

Continuous oxygen concentration monitoring is recommended for patients who require increased oxygen concentration to ensure that the correct oxygen supply is provided and suitable blood oxygen saturations are achieved as prescribed. This device is designed with a connector for pure oxygen. Use of any other gas or mixture of gases will prevent the device from functioning correctly.

- . Keep device away from sources of ignition and/or open flames. Do not smoke near the device.
- . Oxygen concentration and flow of delivered gas should only be determined by a physician or trained medical professional.
- . Oxygen supply shall be standard dry and clean compressed oxygen for medical applications.
- . Shut off the oxygen supply before turning off the device.
- . Before turning on the oxygen supply, run the device to ensure target flow rate can be reached. Only turn on the oxygen supply when the device is running.
- Oxygen leakage may cause serious accidents.
- . If an oxygen leak is discovered, shut off the supply of oxygen and contact authorized personnel for repair. Inspect the oxygen tube and ensure it is free of folds or kinks.




9) Adjusting settings during operation. When the unit is running, "Dew Point Temperature", "Flow Rate", "Oxygen Concentration Alarm Upper Limit", and "Oxygen Concentration Alarm Lower Limit" can be adjusted. The mode (adult/paediatric) cannot be changed.

Step 1:

Rotate the control knob to highlight the lock icon . Press control knob to unlock .


Step 2:

Rotate the control knob to highlight the parameter value (green).



Step 3:

Press the control knob to enter adjustment mode (yellow).



Step 4:

Rotate the control knob to adjust the parameter, then press the control knob to confirm (green). If the control knob is not pressed within 10 seconds after adjustment, the changes made will be discarded, and the system will continue to run with current settings.



Step 5:


Follow the steps 1~3 to change dew point temperature or oxygen concentration alarm limits. If no action has been made for 10 minutes, the device will enter screen saver mode and show the real-time flow rate, dew point temperature and oxygen concentration.

3.3 Humidification Treatment


Connect the patient interface (nasal cannula) to the heated breathing circuit/integrated heated circuit. Read the separate user instructions for the exact patient interface that will be used, including all warnings.

- Warning: Do not bend the breathing circuit or block the patient interface.
- Warning: Note that oxygen concentration above 95% will be displayed as 100%. In this situation, the oxygen reading will flash in red and the device will alarm.
- Warning: Note that oxygen concentration below and equal 25% will be displayed as 21%.


3.4 Settings Menu



When the unit stops, rotate the control knob to highlight the settings icon . Press the control knob to bring up the settings menu.



If there is no operation for 10 seconds, the system will return to the main control interface.



• Time Setting

Select and highlight "Time". Press the control knob to enter the time settings page. Set time Y-M-D-H-M-S.




• Device Information

Select and highlight "Info". Press the control knob to show device information including serial number, model number and software version.





• Choose Language


Select and highlight "Language". Press the control knob to enter the language selection page, and select the desired language.



• Volume Control

Rotate the control knob to highlight the volume  icon, press the control knob to select. Rotate the control knob to adjust the volume setting, then press control knob to confirm.

 Warning: Do not set alarm volume too low when background noise is high!




• History

Select and highlight "History". Press the control knob to enter the history page. Select "Events" or "Alarm History" from this page. Press the control knob to view the corresponding item.



3.5 Alarms


The inspired™ O2FLO is equipped with alarm indicator light, auditory alarm signal and on-screen information signal to warn about interruptions to patient treatment. When an alarm is triggered, the status LED flashes in red or yellow, and an audible alarm will sound. Information describing the alarm will appear on the LCD.

In Events page, up to 100 most recent records of change in settings or operating status can be retrieved.



In Alarm History page, up to 100 most recent alarm records can be retrieved.



 Warning: This alarm system provides the user with the option of adjusting the volume of alarms. If the alarm volume is set too low, the user may not hear it due to high ambient noise level.

• Alarm Alert List


No.	Alarm Information	Alarm Priority	Meaning	Silence Able	Delays	Action
1	Overheat	High	The dew point temperature exceeds 43°C.	Yes	< 5s	Press the pause button to stop the unit running. Wait for a few minutes before restarting the unit. <ul style="list-style-type: none"> •Check the breathing circuit and patient interface for blockage. •Do not change settings or pause too often, especially from high temp. to low temp. and from high flow to low flow.
2	System failure (E#)	Medium	E1 Oxygen concentration transducer failure; E2 Flow sensor failure; E4 Heater plate temp. sensor failure; E5 Cooling fan failure; E6 Turbine fan failure; E8 Heater plate open circuit; E9 Memory failure.	NO	< 10s	<ul style="list-style-type: none"> • Replace a new temp. probe/heater wire/ heated breathing circuit. • Turn off the unit and then restart. If the error still exists, mark down the fault code and contact VM or your supplier.
3	Check chamber	Medium	Using wrong chamber or the chamber is not installed properly.	Yes	< 5s	Use the specified chamber and install the chamber correctly into place.
4	Lack of water	Medium	The chamber water level is lower than the center of prism.	Yes	< 25s	Replace with new water bag with sterile water.
5	Temp. probe uninstalled /failure	Medium	Temperature Probe is not installed properly or has failed.	Yes	< 5s	<ul style="list-style-type: none"> • Install the temperature probe correctly or replace it with a new one. • Replace a new Heated breathing circuit.
6	Check heater wire	Medium	The heater connection cable is not installed properly.	Yes	< 5s	<ul style="list-style-type: none"> • Check if the heater connection cable is connected correctly. • Replace with a new Heated Breathing Circuit. • Replace with a new Heater connection cable.
7	Circuit Failure /data cable uninstalled	Medium	The VUN-300 connection cable is not installed to circuit or the intergrated heated circuit has failed.	Yes	< 5s	<ul style="list-style-type: none"> • Check if the VUN-300 connection cable is connected correctly. • Replace with a new intergrated heated circuit. • Replace with a new VUN-300 Connection cable. • Restart the device. If the problem persists, please contact Vincent medical or your Vincent medical supplier.

No.	Alarm Information	Alarm Priority	Meaning	Silence Able	Delays	Action
8	Incompatible circuit type	Medium	Adult Mode selected when peadiatric integrated heated breathing circuit is connected, or vice versa.	Yes	< 15s	Check if the breathing circuit matches with the selected mode.
9	Check for blockages	Medium	The flow rate is lower than 50% of the target flow.	Yes	< 30s	<ul style="list-style-type: none"> • Check the patient interface or breathing circuit if it is bent, kinked or blocked. • If the inner gas pathway is blocked, contact the professional maintenance personnel. • Unspecified cannula. Please replace with Vincent Medical specified cannula listed on page 11 of this User Manual. • Restart the device. If the problem persists, please contact Vincent medical or your Vincent medical supplier.
10	Oxygen too high	Medium	The oxygen concentration is higher than the upper alarm limit.	Yes	< 20s	Adjust the flow of oxygen from the oxygen supply.
11	Oxygen too low	Medium	Oxygen too low or oxygen cannot reach 90%.	Yes	< 20s	<ul style="list-style-type: none"> • Check if the oxygen gas pathway is blocked or oxygen source is correctly connected. • Adjust the flow of oxygen from the oxygen supply. • Oxygen gas path is blocked. Please repair device.
12	Flow out of range	Medium	The real-time flow is beyond flow specification requirement. The oxygen flow rate is higher than target flow rate.	Yes	< 30s	<ul style="list-style-type: none"> • Check the breathing circuit and the patient interface for blockage. • Adjust the flow of oxygen from the oxygen supply. • Restart the device. If the problem persists, please contact Vincent Medical or your Vincent Medical supplier.

No.	Alarm Information	Alarm Priority	Meaning	Silence Able	Delays	Action
13	Temp. out of range	Medium	After warm up, the chamber outlet dew point temperature is 2°C higher or lower than the target temperature.	Yes	< 30min	<ul style="list-style-type: none"> If the dew point temp. is 2°C below the target temp., decrease the target flow rate. Adjust the ambient temp. by air condition if necessary. If the dew point temp. is 2°C above the target temp., allow a few minutes for the unit to cool down. Restart the device. If the problem persists, please contact Vincent Medical or your Vincent Medical supplier.
14	Temp. probe out	Medium	After warm up, temp. probe out from probe port for more than 30s.	Yes	< 30s	<ul style="list-style-type: none"> Insert the Temp. probe to the sensor hole properly. The ambient temperature is too low . Integrated circuit failure. Please replace with new integrated circuit.
15	SD card uninstalled	Medium	While in operation, the SD card is pulled out.	Yes	< 5s	Insert the SD card properly.

• Alarm mute

When the device alarms, press the control knob to mute the alarm for maximum 120s. Press the control knob again to unmute the alarm.




* This device can record up to 10,000 alarm history records. If necessary, please contact Vincent Medical or your local supplier to obtain the data record software to read the alarm history.


* When the alarm occurs, please follow the action in Alarm Alert List and Troubleshooting to clear the alarm.

• Checking alarm system functionality


The functionality of the alarm system can be checked during self-test when the device is starting up. Please refer to 3.2 Device Setup.

3.6 After Use

After use, press Stand-by button  to shutdown. Ensure that oxygen supply is shut off. If oxygen is flowing or if oxygen concentration is higher than 25%, the window below will appear.



If oxygen concentration is lower than 25%, the device will enter the cooling down mode for blowing and cooling down the heater plate. After cooling down, the device will shutdown; or press and hold standby button for 2 seconds to shutdown.





Disconnect the power plug from the mains power.

If the unit is unplugged when it is running, the “power out” auditory indication will be activated and will sound for more than 120 seconds. If the power cord is plugged back into the power supply within 2 minutes, the unit will resume running normally.

3.7 Data Record Software

The SD card in the SD card slot allows real time usage data of the machine to be recorded. The SD card can be removed by pressing the SD card in. It will then spring out.



Insert the SD card into a computer and install the inspired™ O2FLO Data Reader.exe software. Then run the Vincent Medical Data Record software .

Language
Import data
Export data

Select start date or end date

Edit patient information

Generate a PDF report

The screenshot shows the software interface with several tabs: Basic Info, Usage Time, Detail Charts, Alarm Info, and Work Info. The Patient Info section is highlighted with an 'Edit' button. A 'Generate Report' button is located at the bottom right.

To begin, select "Set" -> "Import Data", then a folder select dialog will pop up. Find the SD Card data folder "VM_VUN". After importing data, the time setting function will be activated. The longest time span is 60 days.


Edit patient information. Name, age and gender are required.

The screenshot shows the 'Edit' dialog box for patient information. The fields include Name, Age, Gender (radio buttons for male and female), Height (with a unit dropdown set to CM), Weight (with a unit dropdown set to KG), Phone No., and Address. The 'ok' and 'cancel' buttons are at the bottom.


After finishing editing, Patient Info entered will be shown.

The screenshot shows the updated patient information in the software interface. The Patient Info section now displays: Name: Nancy, Age: 28, Gender: female, Height: 175, Weight: 70, Phone No.: 13200000000, and Address: China. The 'Generate Report' button is still present at the bottom.

An example of Usage Time Chart is shown below.



Detail Charts include: "Dewpoint Temperature", "Flow" and "Oxygen Concentration".



Alarm Information and work information

Note: A warning will be displayed on the LCD screen when the log becomes full.

The Alarm Information and Work Information interface includes a table with the following columns: No., Date, Set Adult or Paediatric, Flow, Oxygen Concentration, and Dewpoint Temperature. The table lists 17 entries for the date 2017-11-01.

No.	Date	Set Adult or Paediatric	Flow	Oxygen Concentration	Dewpoint Temperature
1	2017-11-01 00:00:00	adult	60	34	37
2	2017-11-01 00:01:00	adult	60	31	37
3	2017-11-01 00:02:00	adult	60	33	37
4	2017-11-01 00:03:00	adult	60	30	37
5	2017-11-01 00:04:00	adult	60	30	37
6	2017-11-01 00:05:00	adult	61	34	37
7	2017-11-01 00:06:00	adult	60	32	37
8	2017-11-01 00:07:00	adult	60	31	37
9	2017-11-01 00:08:00	adult	60	32	37
10	2017-11-01 00:09:00	adult	60	34	37
11	2017-11-01 00:10:00	adult	60	33	37
12	2017-11-01 00:11:00	adult	60	31	37
13	2017-11-01 00:12:00	adult	61	33	37
14	2017-11-01 00:13:00	adult	60	30	37
15	2017-11-01 00:14:00	adult	60	33	37
16	2017-11-01 00:15:00	adult	60	30	37
17	2017-11-01 00:16:00	adult	60	31	37

For the alarm information and work information, user can export data to .csv file. After exporting data, you will get two .csv files: "SN_Alarm.csv" and "SN_WorkData.csv", where "SN" is the 10-digit serial number of the device. To generate report, click "Generate Report" button.

The Patient Information and Generate Report interface shows patient details for Nancy, age 28, female, and a "Generate Report" button.

Patient Info		Patient Info	
Name:	Nancy	Height:	175
Age:	28	Weight:	70
Gender:	female		
Phone No.:	13200000000		
Address:	China		

Generate Report

This will generate a report with a file name like "patient name_phone No.pdf".



4. Product Maintenance




The device should be properly maintained throughout its expected service life. Visual inspection should be conducted on heater plate, temperature probe, power cord and all accessories in use, on a regular monthly basis, by trained technical personnel familiar with the technical performance of the inspired™ O2FLO. If necessary, damaged parts should be replaced. Check that the filter is clean before using the device.

⚠ Warning: Ensure that the device is powered off and disconnected from power supply before any cleaning or maintenance.

⚠ Warning: When the device is not in use, please keep the inlets clean, and store the device in its original packaging.

4.1 Filter Replacement

The air inlet has a filter which needs to be replaced at least once a month. Once discoloration appears, we would recommend the filter to be changed immediately. The filter should be replaced every 30 days for single-patient use, or after every patient. The replacement steps are as follows:



- Rotate counter-clockwise to open the filter cover.
- Remove the filter.
- Clean the filter cover with fresh water.
- Dry the filter cover and filter enclosure with a clean, dry cloth.
- Replace the old filter with a new filter.
- Rotate clockwise to close the filter cover.

** Before installing the filter, ensure that the filter is in good condition without defects.*

4.2 Main Unit Cleaning

Clean main unit, heater connection cable and VUN-300 connection cable with 75% ethanol solvent to wipe the surface of equipment.

⚠ Warning: Do not immerse the main unit in liquid.


⚠ Warning: Do not use corrosive solvents to clean the surface of the main unit.

4.3 Temperature Probe Cleaning and Disinfection


This section only applies to VUN-200 Temperature Probe. Please refer to section 4.2 for the cleaning of Heater Connection Cable and VUN-300 Connection Cable.

a. Rinse the temperature probe in clean water and a mild non-abrasive detergent, and remove any dirt with a soft brush.

b. Flush the temperature probe with water.



⚠ Warning: Do not immerse the probe connector, which connects to the main unit, in any liquid when washing and soaking the temperature probe.




- Dry the temperature probe with a clean, dry cloth.
- Soak for 3 minutes with 75% ethanol.*
- Dry the temperature probe with a clean lint free cloth.
- Pack the temperature probe in clean plastic bag.

*** The 75% ethanol disinfection solution can be replaced by following disinfectants:**

1. Health Essence Disinfecting Effervescent Tablets

2. CIDEX OPA

*** Follow the instructions supplied by the manufacturer of the disinfectant.**



The temperature probe can withstand about 500 cycles of cleaning and disinfection. If the temperature probe fails, the “Temp. probe uninstalled/failure” alarm will be activated when it is connected to the device. In this situation, the temperature probe should be replaced.

4.4 Troubleshooting

Symptom	Possible Cause	Remediation
The LED does not illuminate	Power plug is not plugged into the electrical outlet securely.	Plug the power plug into the electrical outlet correctly.
	Equipment failure.	Contact professional maintenance personnel.
No gas is delivered when the device is running	The gas path is blocked, the circuit is bent or patient interface is blocked.	Check the patient interface and breathing circuit.
	Equipment failure, the blower or the flow sensor is failure.	Contact professional maintenance personnel.
Leakage sound	The check valves are damaged.	Contact professional maintenance personnel.
	The breathing circuit and patient interface are not connected properly.	Check the patient interface and breathing circuit.
	Equipment failure.	Contact professional maintenance personnel.
Oxygen concentration is limited to low levels	The oxygen flow from the oxygen flow meter is lower than the requirement (60LPM).	Use an oxygen flow meter that can deliver oxygen flow up to 60LPM.

4.5 Maintenance Statement

For any repair or change of components, please contact Vincent Medical or your local authorized supplier.

Vincent Medical will provide device related documents to assist service personnel in parts repair.

Installers and operators must follow the instructions of installation, operation, inspection and maintenance. They must be authorized by Vincent Medical.

Inspections and maintenance must be performed in accordance with the recommended schedule. In any of the following circumstances, Vincent Medical will not be responsible for the safety and reliability of the device performance:

- * Modifying, or repairing the machines without any authority from Vincent Medical
- * Unauthorized components are used
- * Electrical power source is not compatible with local regulations
- * Use of device is not in accordance with the instruction manual

It is recommended that users obtain the following information before performing maintenance or repairs:


- * Nature and scope of maintenance or repairs conducted
- * Changes in the scope of maintenance
- * Maintenance date
- * Name of staff or company providing the maintenance service
- * Signature from the Maintenance operator

4.6 Battery Information

One 3V Lithium/manganese dioxide battery is built into the device to power the clock. Only trained personnel should conduct battery replacement. Incorrect replacement would result in a hazard.

4.7 Main Unit Leakage Test

Install a VHC60 chamber to the device. Connect the spike to a water bag. Tightly cover the chamber outlet with cap. Run the machine ensuring no gas is escaping through the chamber outlet. The stable flow rate displayed on screen should be less than 1 L/min for at least 30s. If the flow rate is higher than 1 L/min your device may have a leak in the gas path. Please contact Vincent Medical, your Vincent Medical supplier or authorized personnel for repair.



5. EMC Information



The electromagnetic compatibility (EMC) of the inspired™ O2FLO is designed according to Medical Electrical Equipment Part 1-2 General requirement on safety (IEC60601-1-2), Collateral Standard: Requirements and Tests for Electromagnetic Compatibility and the device complied with the requirements. In certain circumstances, the unit may affect or be affected by nearby equipment due to the effects of electromagnetic interference. If this should happen, move the unit or the location of the unit causing interference, or alternatively consult your homecare provider.


5.1 Declaration - Electromagnetic Emissions

Guidance and Manufacturer's Declaration - Electromagnetic Emissions		
The model VUN-001 is intended for use in the electromagnetic environment specified below. The customer or the user of the model VUN-001 should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The model VUN-001 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The model VUN-001 is suitable for used in domestic establishment and in establishment directly connected to a low voltage power supply network which supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies	

5.2 Guidance & Declaration - Electromagnetic Immunity

Guidance & Declaration - Electromagnetic Immunity			
The model VUN-001 is intended for use in the electromagnetic environment specified below. The customer or the user of the model VUN-001 should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±8kV contact ±2kV, ±4kV, ±8kV, ±15kV air	±8kV contact ±2kV, ±4kV, ±8kV, ±15kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient / burst IEC 61000-4-4	±2kV for power supply lines ±1kV for Input/output lines	±2kV for power supply lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±0.5kV, ±1kV line to line ±0.5kV, ±1kV, ±2kV line to earth	±0.5kV, ±1kV line to line	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11.	< 5% U _T (> 95% dip in U _T) for 0.5 cycle < 5% U _T (> 95% dip in U _T) for 1 cycle 70% U _T (30% dip in U _T) for 25/30 cycles < 5% U _T (> 95% dip in U _T) for 5/6 sec	< 5% U _T (> 95% dip in U _T) for 0.5 cycle < 5% U _T (> 95% dip in U _T) for 1 cycle 70% U _T (30% dip in U _T) for 25/30 cycles < 5% U _T (> 95% dip in U _T) for 5/6 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the model VUN-001 requires continued operation during power mains interruptions, it is recommended that the model VUN-001 be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE: U _T is the a.c. mains voltage prior to application of the test level.			

5.3 Guidance & Declaration - Electromagnetic Immunity Concerning Conducted RF & Radiated RF

Guidance & Declaration - Electromagnetic Immunity			
The model VUN-001 is intended for use in the electromagnetic environment specified below. The customer or the user of the model VUN-001 should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 Vrms 150 kHz to 80 MHz	Portable and mobile RF communications equipment should be used no closer to any part of the model VUN-001, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = [3,5/V1] \times P^{1/2}$ $d = 1,2 \times P^{1/2}$ 80 MHz to 800 MHz $d = 2,3 \times P^{1/2}$ 800 MHz to 2.7 GHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, ^a should be less than the compliance level in each frequency range. ^b Interference may occur in the vicinity of equipment marked with the following symbol: 
	6 Vrms in ISM and amateur radio bands	6 Vrms in ISM and amateur radio bands	
Radiated RF IEC 61000-4-3	10 V/m 80 MHz to 2.7 GHz	10 V/m 80 MHz to 2.7 GHz	NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies. NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.
	385MHz-5785MHz Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communication equipment (Refer to table 9 of IEC 60601-1-2:2014)	385MHz-5785MHz Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communication equipment (Refer to table 9 of IEC 60601-1-2:2014)	
^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the model VUN-001 is used exceeds the applicable RF compliance level above, the model VUN-001 should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the model VUN-001. ^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3V/m.			

5.4 Recommended Separation Distances Between Portable and Mobile RF Communications Equipment and the Model VUN-001



6. Use and Warranty Period



Do not attempt to disassemble the device without authorization. The warranty period is 2 years from the date of purchase. For shelf-life of accessories please refer to the appropriate manual. During the warranty period, the warranty will be void under the following conditions:

- 1) An error caused by operating the unit in any unprescribed conditions or applications.
- 2) Damage caused by using an improper power supply, improper installation or operations not listed in this manual.
- 3) Damage caused by installation, modification (the technical description is included in the user manual) or repair from unauthorized service engineers.
- 4) Damage caused by natural disasters such as fire, earthquake, power surge, lightning, flood, etc.

When the machine develops a fault, please contact your supplier or Vincent Medical for maintenance.

Customer Service Contact Information

Vincent Medical Manufacturing Co., Ltd.

Flat/RM B2, 7/F., Hang Fung Industrial Building, Phase 2, 2G Hok Yuen Street, Hung Hom, Kowloon, Hong Kong

Tel: +852 2186 1010

E-mail: service@inspired-medical.com

Website: www.inspired-medical.com



Obelis s.a

Bd Général Wahis, 53, 1030 Brussels, Belgium

Tel: +32 2 732 59 54

Fax: +32 2 732 60 03

E-Mail: mail@obelis.net



Recommended Separation Distances Between Portable and Mobile RF Communications Equipment and the Model VUN-001			
The model VUN-001 is intended for use in an electromagnetic environment in which radiated RF disturbances is controlled. The customer or the user of the model VUN-001 can help to prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the model VUN-001 is recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter w	Separation distance according to frequency of transmitter m		
	150kHz to 80MHz $d=1,2 \times P^{1/2}$	80MHz to 800MHz $d=1,2 \times P^{1/2}$	800MHz to 2,5GHz $d=2,3 \times P^{1/2}$
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) accordable to the transmitter manufacturer. NOTE 1 At 80 MHz and 800 MHz. the separation distance for the higher frequency range applies. NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			