



HEPA

FSC No.	504-010
Part No.	51050947
Connector(mm)	22F/15M-22M/15F
Use	Disposable
Media Type	Pleated mechanical HAPA filter
NaCl Filter Efficiency	99.97%
Pressure Drop	2.5cmH ₂ O@30 l/min
Compressible Volume	65ml
Weight	48g
Package	100pcs/carton



Low Resistance bacterial / Viral respiratory filter BSF 103-104

FSC No.	510-014	510-014
Part No.	51008900	51009000
Compressible Volume	65ml	65ml
Pressure Drop at 30LPM	1cmH ₂ O	1cmH ₂ O
Bacterial Filtration Efficiency	99.99+%	99.9999+%
Viral Filtration Efficiency	99.99+%	99.99+%
Casing Material	k-resin	k-resin
Weight	27g	27g
Package	100pcs/carton	100pcs/carton



High performance bacterial / Viral respiratory filter BSF 201-204

FSC No.	504-011	504-011	504-011	504-011
Part No.	50421500	50421700	50452200	50453700
Pressure Drop at 30LPM	1.5cmH ₂ O	1.5cmH ₂ O	1.5cmH ₂ O	1.5cmH ₂ O
Bacterial Filtration Efficiency	99.99994+%	99.99994+%	99.99994+%	99.99994+%
Viral Filtration Efficiency	99.997+%	99.997+%	99.9999+%	99.9999+%
Gas Sampling port	yes	no	yes	no
Weight	33g	33g	33g	33g
package	100pcs/carton	100pcs/carton	100pcs/carton	100pcs/carton

HMEF Series



HMEF 1000 Straight

FSC No.	504-011
Part No.	50497100
Connector	22F/15M-22M/15F
Dead Space	35 ml
Tidal Volume Range	150~1000ml
Pressure Drop	0.35kPa@30l/min
Moisture Output	34.5mg/l @ V _t 1000ml
Bacterial Filtration Efficiency	>99.999%
Viral Filtration Efficiency	>99.99%
Weight	30g
Package	100pcs/carton



HMEF Midi Straight

FSC No.	504-011
Part No.	50400216
Connector	22F/15M-22M/15F
Dead Space	22ml
Tidal Volume Range	60~500ml
Pressure Drop	0.35kPa@30l/min
Moisture Output	32mg/l @ V _t 500ml
Bacterial Filtration Efficiency	>99.999%
Viral Filtration Efficiency	>99.99%
Weight	20g
Package	100pcs/carton



HMEF 1000 Angled

FSC No.	504-011
Part No.	50400218
Connector	22F/15M-22M/15F
Dead Space	44ml
Tidal Volume Range	150~1000ml
Pressure Drop	0.35kPa@30l/min
Moisture Output	34.5mg/l @ V _t 1000ml
Bacterial Filtration Efficiency	>99.999%
Viral Filtration Efficiency	>99.99%
Weight	30g
Package	100pcs/carton



HMEF Midi Angled

FSC No.	504-011
Part No.	50400217
Connector	15M-22M/15F
Dead Space	25ml
Tidal Volume Range	60~500ml
Pressure Drop	0.35kPa@30l/min
Moisture Output	32mg/l @ V _t 500ml
Bacterial Filtration Efficiency	>99.999%
Viral Filtration Efficiency	>99.99%
Weight	20g
Package	100pcs/carton

HME Trach



HME Trach Basic

FSC No.	510-019
Part No.	51050901
Connector	15F
Dead Space	14ml
Pressure Drop	0.15 kPa@30l/min
Moisture Output	26mg/l @ V _t =500ml
Weight	15g
Package	200pcs/carton



HME Trach T

FSC No.	510-019
Part No.	51050902
Connector	15F
Dead Space	8ml
Pressure Drop	0.15kPa@30l/min
Moisture Output	26mg/l @ V _t =500ml
Weight	10g
Package	200pcs/carton